
FINNPIPETTE COLOUR

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A Thermo Electron Business

ENGLISH

PRODUCT DESCRIPTION	3
PIPETTE OPERATION	4
PIPETTING TECHNIQUES	4
MAINTENANCE	7
CALIBRATION	8
PACKAGE	11
TROUBLE SHOOTING	11
FIGURES	32
SPARE PARTS	36-37

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PRODUCT DESCRIPTION

The FINNPIPETTE is a continuously adjustable, general purpose micropipette for sampling and dispensing accurate liquid volumes.

It operates on an air displacement principle (i.e. an air interface) and uses detachable, disposable tips.

The adjusted delivery volume is displayed digitally on a readout window in the handle.

FINNPIPETTE covers a volume range of 0,5-5000 μ l with five different models.

Order No.	Volume Range		
4027640	0,5 μ l	to	10 μ l
40270280	5 μ l	to	50 μ l
40270290	20 μ l	to	200 μ l
40270300	100 μ l	to	1000 μ l
4027680	1 ml	to	5 ml

RAW MATERIALS

The FINNPIPETTE is made of mechanically durable and chemically resistant materials.

DESCRIPTION OF TIPS

Finntips are recommended for use with the FINNPIPETTE.

They are made of natural coloured polypropylene, generally regarded as the only contamination free material suitable for tips.

Finntips are also autoclavable (121°C).

Following accessories are useful with handling of pipettes and pipetting of reagents:

Stand	9420290
Shelf hanger	2206040

DIGITAL DISPLAY (Fig.1)

The adjusted delivery volume is clearly indicated by a digital display on the pipette handle.

PIPETTE OPERATION

SETTING THE DELIVERY VOLUME (Fig.2)

1. The delivery volume is set using the operating button on the top of the pipette.

To increase the delivery volume turn the operating button anticlockwise. To decrease turn it clockwise.

2. Ensure that the desired delivery volume clicks into place and that the digits for the volume are completely visible on the handle display.

3. Do not set volumes outside of the pipette's specified volume range! Using excessive force to turn the push button outside the range may jam the mechanism and eventually cause defects.

TIP EJECTION (Fig.3)

Each pipette is fitted with a tip ejector system to eliminate the risk of contamination.

The tip ejector system consists of a tip ejector plunger, spring, spring support and tip ejector bushing.

Direct the pipette towards a suitable waste receptacle and press the tip ejector plunger with your thumb thus disposing the tip safely.

PIPETTING TECHNIQUES

FIGURES 4-7:

A = READY POSITION
B = FIRST STOP
C = SECOND STOP

The operation of the pipette is controlled by the operating button. For the best possible accuracy:

Push and release the button slowly at all times, particularly when working with liquids that have a high viscosity. Never allow the operating button to snap back. Ensure that a clean tip is firmly pushed on to the tip cone of the pipette and there are no foreign particles in the tip itself.

Wet the newly attached tip with the solution being pipetted before any actual pipetting takes place.

This is done by filling and emptying the tip 2-3 times.

Hold the pipette vertically during the liquid intake with the "Grippy" resting on the index finger.

The temperature of the tip and pipette should be equalized to that of the solution.

FORWARD TECHNIQUE (Fig.4)

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button to the first stop.
2. Dip the tip under the surface of the liquid in the reservoir about 1cm deep and slowly release the operating button. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
3. Deliver the liquid by gently depressing the operating button to the first stop. After a delay of about one second, continue to depress the operating button all the way down to the second stop. This action will empty the tip.
4. Release the operating button to the ready position.

If necessary, change the tip and continue with the pipetting.

REVERSE TECHNIQUE (Fig.5)

The reverse technique is suitable for dispensing liquids having high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes.

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.
2. Dip the tip under the surface of the liquid in the reservoir about 1cm deep and slowly release the operating button. This action will fill the tip. Withdraw the tip from the liquid touching it against the edge of the reservoir to remove excess liquid.
3. Deliver the preset volume by gently depressing the operating button to the first stop.

ENGLISH

Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.

4. The remaining liquid is either discarded with the tip or pipetted back into the container.

REPETITIVE TECHNIQUE (Fig.6)

The repetitive technique offers a rapid and simple procedure for repeated delivery of the same volume.

Fill a clean reagent reservoir with the liquid to be dispensed.

1. Depress the operating button all the way down to the second stop.

2. Dip the tip under the surface of the liquid in the reservoir about 1cm deep and slowly release the operating button. This action will fill the tip. Withdraw the tip from the liquid touching against the edge of the reservoir to remove excess liquid.

3. Deliver the preset volume by gently depressing the operating button to the first stop. Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.

4. Continue pipetting by repeating the steps 2 and 3.

PIPETTING WHOLE BLOOD (Fig.7)

(e.g. deproteinization in blood glucose determination)

Use forward technique procedures 1 and 2 to fill the tip with blood.

Wipe the tip carefully with a dry clean tissue.

1. Immerse the tip into the reagent and depress the push button to the first stop, making sure the tip is well below the surface.

2. Release the push button slowly to the ready position.

The tip has been filled with reagent. Do not lift the tip out of the solution.

3. Depress the push button to the first stop and release slowly.

Keep repeating this procedure until the interior wall of the tip is clear.

4. Finally, depress the push button all the way to the second stop to completely empty the tip.

MAINTENANCE

When the FINNPIPETTE is not in use make sure it is safely stored in a vertical position. We recommend FinnpiPETTE stand for this purpose.

SHORT TERM CHECKING

The pipette should be checked at the beginning of each day for dust and dirt on the outside surfaces of the pipette. Particular attention should be paid to the tip cone. No other solvents except 70% ethanol should be used for cleaning the pipette.

LONG TERM MAINTENANCE

0,5-1000 µl (Fig.8)

If the pipette is used daily it should be checked at least twice a year. The servicing procedure starts with disassembly of the pipette.

1. Depress the tip ejector pusher.
2. Insert the tooth of the service tool into the opening at the base of the tip ejector pusher.
3. Pull away the tip ejector shaft and tip ejector pusher.
4. Remove the tip cone by turning it counterclockwise with the service tool.
5. Pull out the piston.
6. Remove the O-ring from the tip cone. In model 5-50 µl the O-ring is located deep inside the tip cone. It can be removed by using the thin end of the piston.
- In model 0,5-10 µl the O-rings cannot be changed. We recommend to regrease the piston only.
7. Clean the piston, the piston spring and the O-ring with a dry napless cloth.
8. Check the cylinder for foreign particles.
9. Grease the cleaned parts with lubricant provided in the package.
10. Reassemble the parts.

1-5 ml (Fig.9)

1. Insert the tooth of the service tool into the opening at the base of the tip ejector pusher.
2. Pull away the tip ejector and tip ejector pusher.
3. Remove the cover.
4. Remove the cylinder by pressing with maintenance pliers.

ENGLISH

5. Pull out the piston.
6. Remove the O-ring from the piston.
7. Clean the piston, the piston spring and the O-ring with a dry napless cloth.
8. Check the cylinder for foreign particles.
9. Grease the cleaned parts with lubricant provided in the package.
10. Reassemble the parts.
11. Check the calibration according to instructions.

CALIBRATION

All Finnpiettes are factory calibrated and adjusted to give the volumes as specified with distilled or deionized water. Normally, the pipettes do not need adjustment, but they are constructed to permit recalibration and adjustment for liquids of different temperature and viscosity.

DEVICE REQUIREMENTS AND TEST CONDITIONS

An analytical balance must be used. The scale graduation value of the balance should be chosen according to the selected test volume of the pipette:

Volume range	readable graduation
under 10 μl	0.001 mg
10-100 μl	0.01 mg
above 100 μl	0.1 mg

Test liquid: Water, distilled or deionized, "grade 3" water conforming ISO 3696. Tests are done in a draft-free room at a constant ($\pm 0.5^{\circ}\text{C}$) temperature of water, pipette and air between 20°C to 25°C . The relative humidity must be above 55%. Especially with volumes under 50 μl the air humidity should be as high as possible to reduce the effect of evaporation loss. Special accessories, such as the evaporation trap, are recommended.

CHECKING THE CALIBRATION

The pipette is checked with the maximum volume (nominal volume) and with the minimum volume or 10% of maximum volume, whichever is higher.

E.g. Finnpiette 0.5-10 μl is tested at 10 μl and 1 μl . A new tip is first pre-wetted 3-5 times and a series of ten pipettings are done with both volumes.

A pipette is always adjusted for delivery (Ex) of the selected volume. Measuring volumes taken from balance is not allowed. If the calculated results are in the limits, the calibration of the pipette is correct.

Procedure:

1. Do 10 pipettings with the minimum volume.
2. Do 10 pipettings with the maximum volume.
3. Calculate the accuracy (A) and precision (cv) of both series.
4. Compare the results to the limits in the Table 1.

If the results are in the limits of Table 1, then the calibration of the pipette is correct. Otherwise the pipette must be adjusted and checked again.

Range	Volume μl	Accuracy		Precision	
		μl	%	s.d.* μl	cv%
0,5-10μl	10	±0.10	±1.0	0.08	0.8
	1	±0.035	±3.5	0.03	3.0
5-50μl	50	±0.30	±0.6	0.15	0.3
	5	±0.15	±3.0	0.13	2.5
20-200μl	200	±1.2	±0.6	0.6	0.3
	20	±0.6	±3.0	0.3	1.5
100-1000μl	1000	±5.0	±0.5	2.0	0.2
	100	±1.5	±1.5	0.6	0.6
1-5ml	5000	±25.0	±0.5	10.0	0.2
	1000	±15.0	±1.5	5.0	0.5

ADJUSTMENT (Fig.10)

Adjustment is done with the service tool.

1. Place the service tool into the openings of the calibration nut at the top of the handle.
2. Turn the service tool clockwise to increase, or counterclockwise to decrease the volume.
3. After adjustment check the calibration according to the instructions above.

FORMULAS FOR CALCULATING RESULTS

Conversion of mass to volume

$$V = (w + e) \times Z$$

V = volume (μl)

w = weight (mg)

e = evaporation loss (mg)

Z = conversion factor for mg/μl
conversion

ENGLISH

Evaporation loss can be significant with low volumes. To determine mass loss, dispense water to the weighing vessel, note the reading and start a stopwatch. See how much the reading decreases during 30 seconds (e.g. 6 mg = 0.2 mg/s). Compare this to the pipetting time from taring to reading. Typically pipetting time might be 10 seconds and the mass loss is 2 mg (10s x 0.2mg/s) in this example. If an evaporation trap or lid on the vessel is used the correction of evaporation is usually unnecessary. The factor Z is for converting the weight of the water to volume at test temperature and pressure. A typical value is 1.0032 µl/mg at 22°C and 95 kPa. See the conversion table on page 39.

Accuracy (systematic error)

Accuracy is the difference between the dispensed volume and the selected volume of a pipette.

$$A = \bar{V} - V_0$$

A = accuracy
 $\bar{V}$ = mean volume
 V_0 = nominal volume

Accuracy can be expressed as a relative value:

$$A\% = 100\% \times A / V_0$$

Precision (random error)

Precision refers to the repeatability of the pipettings. It is expressed as standard deviation (s) or coefficient of variation (cv)

$$S = \sqrt{\frac{\sum_{i=1}^n (V_i - \bar{V})^2}{n-1}}$$

s = standard deviation

$\bar{V}$ = mean volume

n = number of measurements

cv is the relative value of standard deviation.

$$cv = 100\% \times s / \bar{V}$$

CAUTION!

The Finnpiquette is designed to allow easy in-lab service. If you, however, want to send the pipette to us or to our local representative for service, please, enclose a list of any infectious, radioactive or otherwise hazardous materials that have been pipetted.

Also, please, note that the postal authorities in your country may limit the sending of contaminated material by mail.

PACKAGE

The FINNPIPETTE is shipped in a custom made, cardboard package. The package contains the following items:

1. The pipette
2. Service tool
3. Sample of Finntip
4. Tube grease
5. Instruction Manual

TROUBLE SHOOTING

The accompanying table is a guide to possible problems and their solutions.

Defect	Possible reason	Solution
Leakage	Tip incorrectly attached	Attach firmly
	Foreign bodies between tip and tip cone	Clean tip cones attach new tips
	Foreign bodies between the piston, the O-ring and the cylinder	Clean and grease O-ring and cylinder. Use grease
	Insufficient amount of grease on cylinder and O-ring	Grease accordingly
Inaccurate dispensing	O-ring damaged	Change the O-ring
	Incorrect operation	Follow instructions carefully
	Tip incorrectly attached	Attach firmly
Inaccurate dispensing with certain liquids	Calibration altered: caused e.g. by misuse	Recalibrate according to instructions
	Unsuitable calibration. High viscosity liquids may require recalibration.	Recalibrate with the liquids in question.

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Product specifications are subject to change without prior notice. Finnpipette® and Finntip® are registered trademarks of Thermo Labsystems Oy.
Consultez LABSYSTEMS ou votre Distributeur Agréé pour une information détaillée et les tarifs du service L.I.S.A.